



NANTONG FUHUA MEDICAL PACKING CO.,LTD

1NO.A AREA,XINKAI INDUSTRY PARK,XINGFU ROAD,NANTONG DEVELOPMENT ZONE CHINA
TEL: 0086*513*68221299/85171438 FAX:0086*513*85706353

Sterilization Reel Pouch in French Medical Paper

- A: Materials:** Medical grade paper in 60gsm + Medical grade film in 55gsm
- B: Sterilization :** ETO(55 °C) / Steam(121°C / 134°C)
- C: Application:** for packing medical devices, disposables, instruments in hospitals, clinic and laboratory.
- D: Characteristics:**
- Easy to use, can be cut to any length as your needs.
 - With ETO/Steam chemical indicators, more convenient for checking sterilization status.
 - Inner contents can be viewed directly through the transparent film.
- E: Materials Advantages**
- **Excellent Sealing Strength:** Paper fibers are shaped in staggered form, which ensures an excellent heat sealing strength performance, also without peeling directions.
 - **Averaged viscose coated performance** makes good peeling effects, without paper fibers after peeling.
 - **Sufficiently germ-proof effects** depend on it's special construction.
 - **Good Air Permeance** makes a sufficient sterilization effects.
- F: All products** are made in 100,000 purified workshop as ISO13485 standards.

